

Agarose, Low Melting

Catalog Number: 2002
Product of Spain

Applications

- DNA/RNA recovery from gel fragments.
- Tissue and cell culture.
- Viral plaque assays.

Features

- Lower melting and gelling temperatures than standard agarose
- Higher clarity than standard agarose gels.
- Greater sieving capacity.
- Higher gel strength than some of other brands. Easy handling.

Specifications

Moisture	≤ 7 %
Ash	≤ 0.4 %
EEO	≤ 0.12
Sulfate	≤ 0.12 %
For 1.5 % Concentration	
Clarity (NTU)	≤ 4
Gel Strength (g / cm ²)	≥ 500
pH in gel	6.74
Gelling Temperature (°C)	24-28
Melting Temperature (°C)	≤ 65.5
Gel Background	Very Low
DNase / RNase Activity	None detected

LabExpress

---we make your research more efficient

Tel & fax: 734-761-8148
e-mail: info@lab-express.com
<http://www.lab-express.com>